BE UNCOLORED

PROERASERTM

— TATTOO REMOVAL —

PROMOITALIA
WELLNESS RESEARCH

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Proeraser Tattoo Removal is the new laser manufactured and distributed by Promoitalia for removing tattoos and benign pigmented lesions.

Two different wave lengths:

- one for treating dark pigments (black, blue and brown) of professional tattoos, which are removed by using the Nd:Yag (1064nm) laser, that relatively penetrates the skin in depth but in a selective manner
- the other wave length doubled by a KTP (532nm) is used for other colours of professional tattoos (red, orange and some shades of violet) and works on hyper pigmentations such as ephelides, hyperchromatic stains and age spots.

How it works?

Extremely brief pulses of highly concentrated luminous energy are directed on the ink of the tattoo (Chromophore).

When the correct wavelength and a sufficient fluence (energy) are used, the ink of the tattoo is quickly heated and breaks into small fragments, the immune system of the body will see to then remove these particles in shatters, allowing the tattoo to fade.

Safety and effectiveness of tattoo removal by laser

Our laser is efficient for removing undesired tattoos, reducing side effects on patients to minimum.

The best manner to do so is by selecting the correct laser wavelength.

The wavelengths will be absorbed by the ink and not by other skin chromophores.

This selectivity allows the patient to see immediate results.

TECHNICAL SPECIFICATION

Light source	Q-switched ND:YAG solid-state laser
Wavelength	1064nm/532nm
Energy level	400mj-1100mj
Fluency	$\leq 12\text{mj}/\text{cm}^2$
Pulse duration	6ns
Frequency	1,2,3,4,5,6 Hz
Cooling system	Built-in distilled water-circle system
Power requirement	single-phase 220V/10A 5./60 Hz
Dimensions	Portable
Net Weight	12Kg

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Ergonomic probe with
two optical spacers

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Duration and number of laser sessions necessary

For most tattoos, 2 to 5 sessions are necessary in order to remove a tattoo, with sessions from 1 to 3 weeks from each other.

Variations in the number of sessions necessary, are mostly due to specific factors of the tattoo (density, depth, ink type and color), although an immune system compromised in some parts of the body renders the removal of the latter more challenging.

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Q - s w i t c h i n g

The technique denominated q-switching is necessary to produce peak power elevated to ultra-rapid pulses (nanoseconds), indispensable to avoid the formation of scars.


Q-switching implicates the use of a gated electronic mirror and other components in order to allow light to avoid the optical cavity of the laser for only a fraction of a second. Concerning our lasers, these pulses only last a few nanoseconds (billionth of a second), allowing the emission of a very strong pulse, sufficient to heat the ink of the tattoo at different hundredth degrees centigrade for this small fraction of a second.


Without the q-switching technique activated for the emission of laser energy, peak power such to eliminate the ink of the tattoo may not be generated.

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